



BBP 1020

EN / Double Electric 2-Phase Breast Pump ■ User's manual

CZ / Elektrická 2fázová odsávačka mateřského mléka ■ Návod k obsluze

SK / Elektrická 2-fázová odsávačka materského mlieka ■ Návod na obsluhu

HU / Kétfázisú elektromos anyatej leszívó ■ Használati útmutató

PL / Laktator elektryczny dwufazowy ■ Instrukcja obsługi

RU / Двухфазный электрический молокоотсос ■ Руководство по эксплуатации

LT / Dviejų etapų elektrinis pientraukis ■ Naudojimo instrukcija



EN	3
CZ	17
SK	31
HU	45
PL	61
RU	77
LT	91

Table of Contents

Introduction	4
Safety instructions and important notes	5
Description	7
Assembling and disassembling the breast pump	9
Assembling and disassembling feeding bottles	10
Power source and its use	10
Using the breast pump	10
Controlling the pump unit	11
Tips and advice for using the breast pump	12
Storing and heating up breast milk	12
Troubleshooting	13
Cleaning and sterilization	13
Storage	14
Technical specifications	14
Explanation of symbols used on the product, its accessories and packaging	15
Instructions and information about the handling of used packaging	15
Warranty conditions	16

Introduction

- It is generally accepted that breast milk is the natural diet for an infant providing the necessary nutrients for growth and protection against infections.
- The breast pump was developed primarily for the purpose of being able to provide breast milk to the infant even in situations when direct breastfeeding is not possible, e.g. due to sore or cracked breast nipples, or when the mother cannot breastfeed in regular intervals, or in order to help alleviate pain from excessively full breasts from which some breast milk must be expressed.
- The BBP 1020 breast pump with the option of expressing breast milk from both breasts simultaneously is also a great solution for mothers with low breast milk production. Expression from both breasts simultaneously not only improves lactation of breast milk but also saves time.
- The two-phase system of the BBP 1020 breast pump simulates the baby's natural sucking action. In the first stage, the fast suction cycle is emulated, by which the breast nipple is stimulated to induce the lactation of breast milk from the mammary gland. This is followed by the second phase where a deeper slower suction action of the infant is emulated. In this phase, breast milk is expressed from the breast nipple and pumped into a bottle.
- In the stimulation mode, the vacuum intensity can be regulated in 4 levels and in the expression mode in 5 levels. The "Stimulation" button enables the special stimulation mode to be activated at any time during operation, providing the most realistic simulation of the baby's natural sucking rhythm.
- The clearly laid out LCD display shows the operating mode, vacuum intensity, operating time and an indicator signalling the use of the breast pump on the left or right breast or on both breasts simultaneously.
- The silicone massage breast pad is equipped with a massaging function that effectively stimulates the breast to evoke the breast milk lactation process, whilst designed to prevent reverse milk flow when it is pumped out.
- The basic accessories for the breast pump include two feeding bottles and adaptors enabling the use of narrow-neck baby bottles.
- Thanks to its compact dimensions and quiet operation, you can always have the breast pump on hand.


Safety instructions and important notes

- Before first use, thoroughly read the user's manual.
- Only use the breast pump in accordance with the instructions included in this user's manual. Not adhering to these instructions may result in incorrect operation of the breast pump or may damage it. Keep this user's manual for future use.
- This breast pump may be handled by children 8 years of age and older or by persons with physical or mental impairments or by inexperienced persons, if they are properly supervised or have been informed how to handle it in a safe manner and understand the potential dangers. Children must not play with the breast pump and its accessories.
- Only use the breast pump with the original AC adaptor and accessories.
- Before connecting the AC adaptor to a power socket, make sure the voltage stated on its rating label corresponds to the voltage in the power socket.
- For hygiene reasons, the breast pump is intended for repeated use by only a single user.
- Do not use the breast pump outdoors, do not expose it to direct sunlight or put it down on an electric or gas stove.
- The breast pump is not intended to be used in the shower or bath.
- Do not use the breast pump in the event of excessive fatigue or sleepiness or while driving a motor vehicle.
- We do not recommend using the breast pump during pregnancy as this could induce childbirth.
- Do not attempt to take the breast pump off the breast while milk is being expressed. First, turn off the breast pump and only thereafter release the breast from the breast pad.
- In the event that you feel an unpleasant or painful pressure while the milk is being expressed, turn off the breast pump and take the breast out of the breast pad. Perform a gentle massage and heat up the breast in the shower and then repeat the expression process. If problems persist, please consult a lactation consultant or your doctor.

- Do not leave the breast pump without supervision when it is connected to a power socket.
- After finishing using the breast pump, always disconnect it from the power socket.
- Clean and sterilize the breast pump according to the instructions contained in chapter "Cleaning and sterilization". Never submerge the pump unit or AC adaptor in water or any other liquids or sterilize these parts because they do not come into contact with the breast milk or the breast.
- Cleaning may be performed by children if they are older than 8 years of age and under supervision. However, they must not perform sterilization due to the risk of scalding.
- Do not use the breast pump if it is in any way damaged, if its AC adaptor is damaged or if any part that must not be submerged in water fell into water or another liquid.
- Do not attempt to repair the breast pump yourself. Have all repairs or adjustments performed at an authorised service centre. By tampering with the appliance you risk voiding your legal rights arising from unsatisfactory performance or quality warranty.
- We recommend saving the original box, packaging material, receipt and responsibility statement of the vendor or warranty card for the duration of warranty. It is appropriate to move the breast pump in its original box.

Description

Description of accessories




There are two sets of parts A①–⑪ and A⑭–⑲ and a single set of parts A⑫–⑬ included in the original packaging of the breast pump.

Furthermore, the original packaging of the breast pump includes the following replacement accessories: 2 cylinders A②, 2 silicone sealants of cylinder bottom cover A④, 2 valves A⑨, 2 flaps A⑪.

All parts except for the pump unit A⑫ and the AC adaptor A⑬ are made from materials that do not contain BPA (Bisphenol A).

Description of the pump unit



Assembling and disassembling the breast pump

- The packaging of the product was designed to provide the best possible protection to its contents against mechanical damage. Therefore, the individual parts of the product are seated in compartments of the inner packaging that prevent their movement inside the packaging during transport or handling. To remove the individual parts out of the internal packaging, it is necessary to first gently push from underneath the compartment in which the part is located. This will release it from the compartment and it is then easy to take it out.
- Before you start assembling the breast pump, clean and sterilize the individual parts according to the instructions contained in the chapter "Cleaning and sterilization" and thoroughly wash your hands.



Attention:

Never submerge the pump unit A¹² or AC adaptor A¹³ in water or any other liquids or sterilize these parts because they do not come into contact with the breast milk or the breast.

When assembling the breast pump, there must be no water inside the cylinder cover A¹, cylinder A², cylinder bottom cover A³, silicone tubing A⁵ or connector A⁸ and A¹⁰, otherwise suction intensity could be reduced.

- Check that no part of the breast pump is damaged or worn out and that all parts that are used during the assembly are dry. Damaged or worn out parts must be replaced for new ones. To order replacement parts, please contact an authorised service centre.
- Always assemble the breast pump according to the following instructions.
- Insert the breast pad A¹⁵ into the plastic brace, which is a part of the breast pump body A⁶, and push it in as far as possible. Use your fingers to push its edge along the perimeter of the plastic brace to ensure it adheres properly.
- Slide the silicone sealant of cylinder bottom cover A⁴ on to the bottom part of the cylinder bottom cover A³ and push it in so that it slides into the groove along the entire perimeter.
- Seat the cylinder bottom cover A³ on to the top part of the breast pump body A⁶ so that the tooth located in the opening of the top part of the breast pump body points towards the cut-out in the bottom part of the cylinder bottom cover A³. Push the cylinder bottom cover A³ against the breast pump body A⁶ and then turn it all the way in the clockwise direction.
- Insert the cylinder A² with the convex part pointing downwards into the cylinder bottom cover A³ and push it against the cylinder bottom cover A³ along its entire perimeter to ensure that it adheres properly.
- Place the cylinder cover A¹ on to the cylinder bottom cover A³ so that the hole for attaching the connector A¹⁰ points in the opposite direction away from the breast pad A¹⁵ and so that the teeth located along the perimeter of the cylinder cover A¹ fit into the grooves located on the perimeter of the cylinder bottom cover A³. Turn the cylinder cover A¹ clockwise all the way. This will firmly connect the cylinder bottom cover A³ and the cylinder cover A¹.
- Connect the connector A¹⁰ to the cylinder cover A¹ using the silicone tubing A⁵. If the connector A¹⁰ is not connected to the tubing A⁵, first connect these two parts.
- Turn the breast pump body A⁶ so that the cylinder cover A¹ points downwards and insert the silicone sealant A⁷ into the breast pump body A⁶. Push it down with your fingers until it slides into the groove along the entire perimeter. There is a tab located on the inner perimeter of the silicone sealant A⁷ that makes it easier to take the sealant off the breast pump body A⁶. When the silicone sealant A⁷ is installed, the tab will be located just above the edge of the groove.
- Connect the flap A¹¹ to the valve A⁹ so that the semicircular-shaped stub on the underside of the flap A¹¹ fits into the hole of the same shape in the valve A⁹. To ensure that the flap A¹¹ properly adheres to the valve A⁹, push it down with your thumb.
- Fasten the valve A⁹ with the flap A¹¹ to the underside of the breast pump body A⁶.
- Turn the breast pump body A⁶ so that the cylinder cover A¹ points upwards and screw the bottle A²¹ into the bottom part of the breast pump body A⁶. Place the bottle A²¹ into the stand A²².
- Connect the free end of the silicone tubing A⁵ with the connector A⁸ to the left or right socket B⁷ or B¹⁴ in the pump unit A¹². If the connector A⁸ is not connected to the tubing A⁵, first connect these two parts.
- If you wish to express milk from both breasts simultaneously, assemble the second breast pump in the same way and connect it to the second available socket on the pump unit A¹².
- Now the breast pump is assembled and ready for breast milk expression.
- To disassemble the breast pump proceed in the reverse sequence.

Assembling and disassembling feeding bottles

- The instruction below describe how to assemble a feeding bottle and how to close it if you wish to store breast milk in it.
- Before you start assembling the feeding bottle, clean and sterilize the individual parts according to the instructions contained in the chapter "Cleaning and sterilization" and thoroughly wash your hands.
- Assemble the feeding bottle by placing the nipple A17 into the bottle opening cap A16 and then screw the bottle opening cap A16 together with the nipple A17 on to the feeding bottle A21. If you wish to use a narrow-neck feeding bottle, connect the adaptor A19 to the bottle opening cap A16 and nipple A17 from underneath and then screw this assembly on to the narrow-neck bottle. If you do not use the feeding bottle immediately after assembling it, we recommend that you place the protective nipple cap A16 on to the bottle opening cap A16 and nipple A17, which will protect the nipple A17 against dust.
- If you wish to store the milk in the bottle, close it according to the following instructions: Turn the bottle opening cap A16 so that the opening with the thread for screwing on the feeding bottle A21 points downwards. Place the bottle cover A20 with the convex side pointing downwards in the centre of the bottle opening cap A16 and push it down thoroughly along its perimeter until an audible click is heard. Turn the bottle cover A20 so that the bottle opening cap A16 is pointing upwards and screw it on to the neck of the feeding bottle A21. Now the feeding bottle A21 is closed. You can check for leaks by turning the feeding bottle A21 with the bottle cover A20 upside down.

Note:

Store and heat the milk according to the instructions contained in chapter Storing and heating up breast milk.

- For disassembly of the feeding bottle proceed in the reverse sequence.

Power source and its use

- Only use AC adaptor A13 model number YHA-24W-05FEU060300 (output 6V $\overline{-}$, 3A) that is supplied with this breast pump.
- Connect the AC adaptor A13 to socket B1 in the pump unit A12. Insert the plug of the AC adaptor A13 into an el. power socket.

Using the breast pump

- Prior to expressing breast milk, thoroughly wash your hands and wipe the breast(s) with a wash cloth (ideally a disposable type) dipped in clean warm water.
- Seat yourself comfortably in a chair or armchair and rest your back against a backrest. It is not appropriate to express breast milk when bending over. Take the breast pump out of the bottle stand A22. Place the breast pad A15 against the breast in such a way that the nipple is located in the middle of the breast pad A15. Use your thumb and index finger to hold the breast pad A15 against the breast and use the palm of the same hand to support the breast. The breast pad A15 must perfectly adhere to the breast. If you wish to express milk from both breasts simultaneously, do the same with the other breast.
- Make sure that the pump unit A12 is connected to a mains socket and start it by pressing button B16. The functions of the control buttons are further described in chapter Controlling the pump unit.
- When you have finished, turn off the pump unit A12 and take the breast out of the breast pad A15. Do not attempt to take the breast out of the breast pad A15 while the breast pump is running. Disconnect the pump unit A12 from the mains power socket.
- Place the feeding bottle A21 with the expressed milk in the bottle stand A22 to prevent it from being tipped over. The screw it out of the breast pump body A6 and screw on the bottle cap A19 with the nipple A17 if you wish to feed the expressed milk immediately, or close it according to the instructions contained in the chapter Assembling and disassembling feeding bottles.

Note:

Do not fill the feeding bottle A⑫ with more than 140 ml of milk so as to prevent spilling it.

- Clean all used parts according to the instructions in chapter Cleaning and sterilization. Clean the feeding bottle A⑫ with the bottle cap A⑬ and nipple A⑭ or bottle cover A⑮ later, only after using up the contents.

Controlling the pump unit

■ Button for turning the pump unit on and off

- To turn on the pump unit A⑫, press button B⑯.
- The breast pump will start automatically in the stimulation mode. The suction intensity is preset to the lowest level 1 and only the left socket B⑦ will be running for stimulation of the left breast. The display will show indicator B④, the operating time B⑤, the corresponding number of segments on the suction intensity scale B⑥ and the stimulation mode indicator B⑧. After 45 seconds of operation, the suction intensity automatically increases by one level and after every 45 seconds of operation it goes up by another level until it reaches suction intensity level 4. The stimulation phase ends automatically after 3 minutes.

Note:

If milk starts flowing out prior to the stimulation phase ending automatically, manually switch the breast pump to the expression mode by pressing button B⑰. The use of button B⑰ is described in more detail in chapter Button for increasing and decreasing suction intensity.

- Upon completion of the stimulation phase, the breast pump automatically switches to the expression mode. In the expression mode, the suction intensity is preset to level 5. The display will show indicator B④, the operating time B⑤, the corresponding number of segments on the suction intensity scale B⑥ and the expression mode indicator B⑨. The breast pump will turn itself off automatically after 30 minutes of operation.
- To turn the breast pump off at any time during operation, press button B⑱.

■ Button for switching between the left and right socket

- Button B⑲ switches between the left and right socket B⑦ and B⑭ for stimulation of the left or right breast or expression from it. The display B② shows either indicator B④ or indicator B⑩.
- Whenever button B⑲ is pressed, the display temporarily shows information about the set suction intensity level in the format of "L + number".

■ Button for starting operation of both sockets

- Button B⑲ serves to put both sockets B⑦ and B⑭ into operation. Pressing button B⑲ once, activates the mode, where sockets B⑦ and B⑭ are switched on in an alternating fashion. Pressing button B⑲ again, activates the mode, where sockets B⑦ and B⑭ are running simultaneously. The alternating operation of sockets B⑦ and B⑭ is signalled by indicators B④ and B⑫ being lit in an alternating fashion. When sockets B⑦ and B⑭ are running simultaneously, both indicators B④ and B⑫ are lit simultaneously.
- Whenever button B⑲ is pressed, the display temporarily shows information about the set suction intensity level in the format of "L + number".

■ Button for increasing and decreasing suction intensity

- Button B⑰ serves to increase suction intensity. Button B⑳ serves to decrease it. Each time button B⑰ or B⑳ is pressed, the suction intensity will be increased or decreased by one level. When the suction intensity level increases, the speed of stimulation/expression automatically decreases, and vice versa.
- In the stimulation mode, it is possible to regulate suction intensity in the range from 1 to 4. Indicator B⑧ is lit on the display. After setting suction intensity level 5, the breast pump automatically switches to the expression mode. This is signalled by the indicator B⑨ being lit. Switch to the expression mode when milk starts flowing into the feeding bottle A⑫. The suction intensity setting in the expression mode can be regulated in the range from 5 to 9.

- If, after starting the pump unit A②, you set the suction intensity in the stimulation mode manually, the breast pump will continue running at the set suction intensity until you change the settings.
 - Whenever button B⑪ or B⑬ is pressed, the display temporarily shows information about the set suction intensity level in the format of "L + number".
- **Button for setting the stimulation mode**
- When button B⑬ is pressed, the breast pump will switch to a stimulation mode, in which it will continuously repeat 10 cycles of stimulation at a suction intensity level of 3 followed by 1 cycle of expression at a suction intensity level of 7. This mode most realistically imitates the infant's natural sucking action in order to achieve the most effective stimulation.
- Note:**
After 3 minutes, if you have not pressed button B⑪ or B⑬ in the interim, the display B② will show indicator B⑮ despite the breast pump still running in stimulation mode.
- Buttons B⑪ and B⑬ can be used to change the preset suction intensity level 3 in the range from 1 to 5. 10 cycles of stimulation/expression will then continuously repeat at the selected suction intensity with 1 suction cycle at a preset suction intensity of 7.
 - The stimulation mode will continue until you switch the breast pump to a standard expression mode. To switch the breast pump to the standard expression mode, it is necessary to first use button B⑪ to increase suction intensity to level 5 and then to press button B③. This will cancel the stimulation mode setting.

Tips and advice for using the breast pump

- If you feel an unpleasant pressure when a high suction intensity level is set, reduce the suction intensity.
- If you are not successful expressing any milk within the time interval of approx. 5 minutes, stop the breast pump. First try relaxing for a while and warm the breasts up in a warm shower. Then try again to express the milk or postpone it to a later time when more milk has been created.
- Usually, 10 minutes should be sufficient to express the milk. This time is only indicative and may vary from person to person.
- If you wish to increase milk production, express from both breasts simultaneously for approximately 15 minutes.

Storing and heating up breast milk

- Only breast milk collected using a clean and sterilized breast pump is suitable for storage. Store the milk in sterilized containers designed for storing breast milk or in special disposable breast milk storage bags.
- To prevent the milk from leaking or absorbing other odours from an area in which it is stored, thoroughly close containers or bags before placing them in the refrigerator or freezer.
- Breast milk can be stored for up to 48 hours in a refrigerator (except for the door area) and up to 3 months in a freezer at a temperature of -18°C .
- If you express the milk gradually and wish to mix the individual batches together into a single container and freeze them, do so no later than within 24 hours. It is only possible to mix batches that are cooled to the same temperature.
- Mark the containers or bags that you place in the freezer with a date.

Note:

A close feeding bottle A⑳ can be stored in a refrigerator or freezer.

- Before this milk is fed to an infant, it should be heated up to a temperature of approx. 37°C . Frozen milk can be left overnight to defrost in a refrigerator. Defrosting in a refrigerator must not, however, take longer than 24 hours. Heat the milk gently in a baby bottle warmer or in a container with warm water. Do not use boiling water as this could destroy the nutrients through high temperatures. Microwave ovens are not particularly suitable for heating breast milk, namely because of uneven heating.

- All parts of the feeding bottle from which the milk is fed to the infant must be sterilized before they come into contact with the breast milk. Before feeding, always first shake the feeding bottle and check the temperature of the milk. Dispose of any unused milk that has already been defrosted or heated up.

Troubleshooting

- In this chapter you will find possible solutions to problems that you may encounter when using the breast pump. If you are unable to remedy the problem according to the following instructions, please contact an authorised service centre.
- **The breast pump does not start when button B⑩ is pressed.**
 - Check that the AC adaptor connector A⑬ is inserted in socket B① on the pump unit A⑫ and that the AC adaptor plug A⑭ is inserted in a working mains power socket.
- **Low or no suction intensity**
 - Check whether the breast pump is correctly assembled and whether the individual parts are properly attached to one another. The connections between the individual parts of the breast pump must not be loose.
 - Check that the breast pad A⑮ has properly adhered to the breast. Do not hold the breast pump by the breast pump body A⑥. While expressing milk, hold the breast pad A⑮ on the breast with your thumb and index finger and support the breast with your palm.
 - Check that the silicone tubing A⑤ is not twisted and that there is no water inside it. During assembly and prior to starting expressing milk, all components of the breast pump must be dry and the silicone tubing A⑤ must not be bent.

Cleaning and sterilization

Cleaning the breast pump

- Make sure that the pump unit A⑫ is turned off and disconnected from the mains power socket. Disassemble the breast pump into its individual parts, see chapter Assembling and disassembling the breast pump

Note:

For the purposes of cleaning, it is not necessary to remove flap A⑪ from the valve A⑨ or the silicone sealant of cylinder bottom cover A④ from the cylinder bottom cover A③.

- **Before first use and after every use**
 - Thoroughly wash the cylinder A②, the breast pump body A⑥, the silicone sealant A⑦, the hygiene cover A⑱, the breast pad A⑮, the nipple cap A⑩, the bottle opening cap A⑰, the bottle adaptor A⑲, the bottle cover A⑳, the feeding bottle A㉑ and the bottle stand A㉒ in warm drinking water with a gentle kitchen detergent. Then thoroughly rinse all these parts under running drinking water and allow them to dry naturally. The above mentioned parts may also be washed in a dishwasher (in the top rack or in the cutlery compartment).

Note:

It is sufficient to wash only the parts that were used after every use. If, for example, you do not use the nipple cap A⑩, the bottle opening cap A⑰ or the bottle adaptor A⑲ immediately, it is sufficient to clean these parts directly before their planned use.

- Carefully wash the valve A⑨ with the flap A⑪ in warm drinking water with gentle kitchen detergent. Then thoroughly rinse under running drinking water and allow to dry naturally.
- Wipe the cylinder cover A①, the cylinder bottom cover A③ and the tubing A⑤ with connectors A⑧ and A⑱ with a clean lightly dampened cloth and allow these parts to dry naturally. These parts do not come into contact with breast milk or with the breast and therefore need not be washed under running water. If you do, however, wash these parts, it will be necessary to wait for the water that entered inside these parts to dry out.

- Never submerge the pump unit A⑫ or the AC adaptor A⑬ in water or any other liquid. If it becomes dirty, you can wipe the surface of the pump unit A⑫ with a clean lightly dampened cloth and then wipe dry.

Sterilizing the breast pump

■ Before first use and every subsequent use:

- Sterilize the cylinder A②, the breast pump body A⑥, the silicone sealant A⑦, the valve A⑨ with the flap A⑪, the hygiene cover A⑭, the breast pad A⑮, the nipple cap A⑯, the bottle opening cap A⑰, the bottle adaptor A⑱, the bottle cover A⑲, and the feeding bottle A⑳ in a steam sterilizer or boil in a pot for approximately 5 minutes. When sterilising by boiling, we recommend using distilled water to prevent undesirable white deposits from forming on the individual parts. Leave the sterilized parts to dry naturally on a clean wiping cloth or clean tray.



Attention:

All parts must be clean prior to sterilization. Cleaning instructions are provided in the preceding subchapter Cleaning the breast pump.

- The bottle stand A㉑, the cylinder cover A①, the cylinder bottom cover A③ and tubing A⑤ with the connectors A⑧ and A⑩ do not need to be sterilized. These parts do not come into contact with breast milk or the breast. If, however, you do sterilize these parts, then it will be necessary to wait for all the water that entered into the cylinder cover A①, cylinder bottom cover A③, tubing A⑤ or connectors A⑧ and A⑩ to dry out completely.
- Do not sterilize the pump unit A⑫ and the AC adaptor A⑬. These parts do not come into contact with breast milk or the breast and sterilization would damage them.

Storage

- Place the cleaned breast pump into, for example, a clean box or wrap it in a clean wiping cloth to protect it from dust. The breast pump can be stored assembled or disassembled. To protect it against dust, you can fit the hygiene cover A⑭ on to the breast pad A⑮.

Technical specifications

Model number	BBP 1020
Voltage and power input	6 V $\overline{\text{---}}$, 1,2 A
Power source	AC adaptor, YHA-24W-05FEU060300 input: 100–240 V $\overline{\text{---}}$, 50/60 Hz, 0,6 A, output: 6 V $\overline{\text{---}}$, 3 A
Protection class against injury by electrical shock	AC adaptor – protection class II
Noise level	60 dB(A)
Suction intensity	Max. –0.02 MPa (stimulation mode) Max. –0.04 MPa (expression mode)
Stimulation/expression speed	96 to 139 cycles/minute (stimulation mode) 38 to 79 cycles/minute (expression mode)
Automatic shut off timer	30 minutes
Operational conditions	Ambient temperature: 5 to 40 °C Relative humidity: \leq 65 %

Explanation of symbols used on the product, its accessories and packaging



This symbol indicates "Important notice".



This symbol indicates that prior to using the product, the user should read the user's manual.



This symbol indicates that this appliance meets all the basic requirements of EU directives related to it.



This symbol indicates that the product has protection class II.



This symbol indicates the polarity of the power cord connectors.



This symbol indicates that it is a product that may be used indoors.



This symbol serves to designate the production lot number.



This symbol indicates that a fee has been paid to an organisation performing recollection and utilisation of the packaging waste.



This symbol indicates the material from which the packaging is made.



This symbol on products or original documents means that used electric or electronic products must not be added to ordinary municipal waste. For proper disposal, renewal and recycling hand over these appliances to determined collection points. Alternatively, in some European Union states or other European countries you may return your appliances to the local retailer when buying an equivalent new appliance. Correct disposal of this product helps save valuable natural resources and prevents potential negative effects on the environment and human health, which could result from improper waste disposal. Ask your local authorities or collection facility for more details. In accordance with national regulations penalties may be imposed for the incorrect disposal of this type of waste.

For business entities in European Union states

If you want to dispose of electric or electronic appliances, ask your retailer or supplier for the necessary information.

Disposal in other countries outside the European Union.

This symbol is valid in the European Union. If you wish to dispose of this product, request the necessary information about the correct disposal method from the local council or from your retailer.

Instructions and information about the handling of used packaging

Dispose of used packaging material at a site designated for waste in your municipality.

Changes to text and technical parameters are reserved.

